


**SEDANA**MEDICAL  
the AnaConDa® technology people



AnaConDa – Anaesthetic Conserving Device

# Inhalation Sedation in Intensive Care Units

Handbook for clinical use



## **Sedana Medical AB Contact Information**

### **Sedana Medical Sweden**

Sedana Medical AB  
Kungsgatan 62  
SE-753 18 Uppsala  
Phone, Switchboard +46 (0)18 18 66 27  
Fax +46 (0)18 10 15 80  
Email [info@sedanamedical.com](mailto:info@sedanamedical.com)

### **Sedana Medical Germany**

Sedana Medical AB  
Am Hofanger 8  
D-82538 Geretsried - Gelting  
Phone +49 (0)8 171 911 341  
Fax +49 (0)8 171 910 072  
Email [robert.vomdorp@sedanamedical.com](mailto:robert.vomdorp@sedanamedical.com)

### **Important user information**

This manual is a technical description of the Sedana Medical AnaConDa<sup>®</sup> system. It has been written to complement the instructions for use and includes detailed descriptions of set-up, medical application and system description. Be sure to carefully read the instructions for use before operating the AnaConDa<sup>®</sup>. Please direct any remarks or comments about this publication to Sedana Medical AB.

### **Trademark**

AnaConDa<sup>®</sup> is a trademark of Sedana Medical AB.

# Table of contents

	Section
<b>1. Introduction</b>	<b>2</b>
Analgesedation of intensive care patients	2
General information about volatile anaesthetics	4
<b>2. AnaConDa – Brief Description</b>	<b>6</b>
<b>3. Technical Requirements</b>	<b>8</b>
Anaesthetic gas scavenging	8
Connecting scavenging systems	11
Anaesthetic gas monitor	11
Syringe pump	12
AnaConDa material/equipment	12
<b>4. Practise Recommendations</b>	<b>13</b>
Filling and positioning of the syringe	13
Start-up and set-up of AnaConDa	14
Brief instructions	15
Dosage calculation for end-tidal concentration (Fet)	16
Practical handling	17
Syringe changes	17
AnaConDa system change	17
Suctioning	17
Using nebulisers	17
Disconnection / connection	18
Patient transport	18
Temporary change in concentration	18
<b>5. Warnings, Notices, Tips and Tricks</b>	<b>19</b>
<b>6. Literature</b>	<b>22</b>
AnaConDa in intensive care	22
Volatile anaesthetics in intensive care units	23
Sedation in general	23
<b>7. Sedation Scoring Systems</b>	<b>24</b>
Ramsey Sedation Score	24
Richmond Agitation Sedation Score (RASS)	25
<b>8. Composition of Materials, Product Specifications, Ordering Information</b>	<b>26</b>
<b>9. Contributing Authors</b>	<b>28</b>

# 1. Introduction

## Analgesedation of intensive care patients

Concepts of analgesedation are quite complex due to the various individual goals of therapy (analgesia, sedation, anxiolysis) and nowadays belong to the standard repertoire of intensive care. The pharmacological development of rapid acting, thus well controllable substances on one hand, as well as knowledge of the connection between analgesedation and outcome based on scientific studies have resulted in a marked change in the premises of analgesedation. While in the past the main focus was primarily on deep sedation, it is now on analgesia, which may be supplemented with sedation adapted to the specific needs (target score RAMSAY -2 to -3). This change in the paradigm is supported by clinical studies in which 70% of patients specified pain as the most unpleasant recollection, despite medical staff having assessed the pain therapy as adequate. This fact concurrently reveals the requirement for adequate monitoring.

In the clinical situation, however, with a sophisticated approach there are still other medical reasons warranting deep analgesedation as part of a



therapy strategy (severe septicaemia, acute phase polytrauma, burn patients, burst suppression for neuroprotection and ICP reduction, etc.). In these cases it is reasonable to routinely check the level of sedation, for example through clinical sedation scores, machine diagnosis (Bispectral-Index, EEG) or daily sedation pauses.

Considering this and the many undesirable side effects during analgo-sedation in clinical routine – such as developing an increasing tolerance, insufficient sedation quality when combining medications (polypragmatism), gastrointestinal immotility, developing withdrawal symptoms, cognitive deficits – make the search for the ideal sedation concept far from over.

**The ideal analgo-sedation method involves:**

1. Good controllability, with quick onset of action and rapid elimination when stopped
2. Predictable effect
3. Organ-independent elimination
4. No accumulation of active metabolites, no enzyme induction
5. Little or no impairment of organ function, particularly the cardio-pulmonary system and the gastrointestinal tract
6. Lowest possible interaction with other medications
7. (No release of histamines)
8. Non addiction with long-term use

Based on the progression of the disease a distinction is made between intensive care stabilization phase and weaning phase. During the acute phase the focus is typically on stabilising the patient, active patient participation is not always necessary. The weaning phase, however, requires active patient participation to wean the patient from the respirator or for physical therapy. Here the focus shifts to analgesia and psycho-vegetative screening (anxiolysis). But in addition behaviour disorders often emerge which can sometimes be traced back to earlier analgo-sedations such as post-operative delirium, excitation. Thus measures should already be taken at the start of analgo-sedation to considerably reduce the emergence of delirium.

## General information about volatile anaesthetics

The use of inhalation anaesthetics is subject to the specific physical principles of absorption, dispersion and elimination. In general the absorption of inhalation anaesthetics implies its gaseous physical condition, in which at room temperature the substances may at first be in liquid form, turning into a gaseous (volatile) state during application. In general anaesthesia, liquid inhalation anaesthetics evaporate through the use of so-called anaesthetic gas evaporator (vapour). The AnaConDa System does not need additional technical requirements to function. It merely utilises the physical properties of the anaesthetic gases Isoflurane and Sevoflurane. Liquid inhalation anaesthetics not only vaporize at their boiling point but partially already at room temperature. In their gaseous stage all inhalation anaesthetics have a specific gas pressure which is substance specific as well as temperature sensitive and significantly impacts absorption in the organism. The higher the gas pressure, the faster the anaesthetic can reach a concentration exchange between alveolar air and blood.

In general the absorption of inhaled anaesthetics depends on:

- Diffusion parameters (alveolar surface, alveolar diffusion distance, etc.)
- Blood solubility
- Cardiac output
- Alveolar pulmonary venous partial pressure difference
- Alveolar concentration

In addition to the density or the molecular weight of the anaesthetic gas the alveolar surface and the diffusion distance (alveolar cell – basal membrane – endothelial cell) are significant clinical determinants for transporting the gas between alveolar air and capillary blood, thus the efficacy of volatile anaesthetics. Under normal conditions the alveolar surface with an area between 55-100 m<sup>2</sup> represents the diffusion area. A reduction of this exchange area, e.g. through atelectasis or emphysema, results in a direct proportional reduction of the diffusion area. The diffusion distance, in contrast, is approx. 0.25-0.6 µm in healthy people and affects the diffusion capacity with the reciprocal of the square of the distance. In clinical applications, doubling the diffusion distance e.g. through a pulmonary oedema, thus reduces the diffusion rate by a factor of 4.

The dispersion of volatile anaesthetics in the human organism is determined by

- Tissue solubility
- Perfusion in the individual organs and tissue

A major portion of anaesthetic gases is eliminated pulmonically. The metabolised portion amounts to 0.2% with Isoflurane and with Sevoflurane this portion is between 3 and 5%. Upon terminating the anaesthetic gas supply a partial pressure gradient quickly develops between the pulmonary vascular bed and the alveolar space, which reversely leads to diffusing the anaesthetic gas from the bloodstream into the alveolar. The flushing of anaesthetic gases is ultimately proportionally impacted by the level of ventilation.

The intensity of a volatile anaesthetic is characterised by the minimal alveolar concentration (MAC). Among other items it also depends on the age, body temperature as well as other sedating and analgesic medications. MAC50 for example is defined as the alveolar concentration of an inhalation anaesthetic at which 50% of all patients no longer respond to incision ("surgical MAC"), MAC95 the concentration at which 95% of patients no longer respond with pain or the MACawake concentration so often aimed at during inhaled sedation at which 50% of patients open their eyes.

## 2. AnaConDa - Brief Description

AnaConDa, an acronym for "Anaesthetic Conserving Device", is a medical device which allows the administration of volatile anaesthetics. It is used along with conventional respirators – without breathing circuit, breathing lime and anaesthesia evaporator – and connected between the Y-piece and endotracheal tube, just as a bacteria/virus filter. It further requires a syringe pump, anaesthetic gas monitor and anaesthetic gas elimination.

The core of the AnaConDa contains an evaporator and a reflector. Altogether the system consists of the following components.





fig. 1:


1. AnaConDa syringe with screw coupling and screw cap
2. Agent supply line with valve and screw coupling
3. Evaporator
4. Gas measurement port
5. Bacteria and virus filter
6. Anaesthesia reflector and humidification filter



The evaporator is a white, hollow bar made from porous plastic. The Ana-ConDa syringe is filled with liquid Isoflurane or Sevoflurane and a syringe pump moves it through the agent line into the evaporator where it leaks onto the surface through pores, evaporating instantly. The anaesthesia reflector consists of activated carbon fibres woven into the cotton-like humidification filter. Anaesthesia molecules attach to the surface of these carbon fibres during expiration, releasing them with the next inspiration. In the process 90% of the anaesthesia molecules are retained in the reflector and re-supplied with the next inspiration. A maximum of 10% of the anaesthesia molecules pass through the reflector and are exhaled via the respirator's expiration valve. This lost amount is continuously replaced by the infused liquid anaesthetics.

### Functional description

Maximum 10% is passing through the Anaesthesia reflector



### 3. Technical Requirements

#### **Anaesthetic gas elimination**

Sedana Medical recommends the use of anaesthetic gas elimination whilst the AnaConDa system is being used.

The AnaConDa (reflector) efficiency is sufficiently high to keep ambient air contamination at approx. 0.5 ppm (ppm = parts per million) under standard intensive care unit conditions (ambient air exchange rate of 8x/h or higher), even without the use of anaesthetic gas elimination.

#### **Anaesthetic gas scavenging**

If an intensive care unit has an active anaesthetic gas scavenging system or vacuum system, it may be utilized for the AnaConDa system. Here the expiratory valve of the intensive respirator is connected to the active gas scavenging system or vacuum via a tube. The use of additional equipment may be required to ensure the flow sensor system on the intensive care respirator is not affected by the negative pressure generated. This accessory is available from the companies Dräger Medical (Lübeck) and Maquet (Sweden).



*EVAC 180 by Maquet*



*AGS by Dräger Medical*



Two mounts are available for the Contrafluran anaesthetic gas filter:



*Sensor mount  
Sensofluran  
Item no. ZE0000052*



*Standard mount  
Item no. ZE0000051*

### **Sensor mount Sensofluran**

This sensor mount includes an ionisation measurement chamber. It measures the fill level of the Contrafluran anaesthetic gas filter and accurately indicates it with LEDs.

### **Standard mount**

Per the manufacturer (ZeoSys, Berlin), the Contrafluran filter has an intake volume of 500 ml liquid anaesthetic (Isoflurane, Sevoflurane). On the standard mount a mark is made on the filter after each injection (50 ml) to control the fill level (10 injections, i.e. 10 marks amount to 500 ml). The filter is then to be replaced.

## Connecting the anaesthetic gas filter to the intensive care respirator



Accessory set  
Item no. 26072

- 22 mm flex tube (Intensive care respirator and anaesthetic gas filter)
- 8 mm O<sub>2</sub> connection tube or 8 mm silicone tube (gas bypass filter and anaesthetic gas filter)
- Adapter 33w/22m connector (Hamilton Galileo, Maquet Servo I, Siemens Servo)
- Intensive care respirator adapter 22m/22w with side outlet to simultaneously feed expiratory air as well as gas monitor expiration air into the anaesthetic gas filter

## Anaesthetic gas monitor

An Anaesthetic gas monitor should be used whilst administering volatile anaesthetics to measure the end-tidal gas concentration (F<sub>et</sub>), providing information on the level of sedation. Two options (main or side-stream flow measurement) are available for measuring the volatile anaesthetic administered (Isoflurane or Sevoflurane and CO<sub>2</sub>).

Be sure to use an anaesthetic gas monitor with moisture separator and/or Nafion tubing (semi-permeable membrane) for bypass flow measurements, as the measuring tube feeds a large amount of moisture to the anaesthetic gas monitor.

Monitors with side-stream flow measurement also require CO<sub>2</sub> or gas measurement tube (possibly Nafion tubing).



Vamos gas monitor  
Vamos by Dräger Medical  
side-stream flow monitor



Anastasia gas monitor  
by Acutronic Medical  
main flow monitor

## Syringe pump


The syringe pump must be set to syringe type

- Becton Dickinson Plastipak 50 ml
- Sherwood Monoject 50 ml

and a maximum cut-off pressure set.

Some syringe pumps require separate settings to be made, but in general any syringe pump may be used.

## AnaConDa material/equipment



AnaConDa System with syringe Item no. 26000



AnaConDa syringe Item no. 26022



Fill adapter  
Isoflurane / Sevoflurane  
Standard thread  
Item no. 26064



Fill adapter Sevoflurane  
Quickfill by Abbott  
Item no. 26042

## Initialisation / AnaConDa set-up

1. Switch on the prepared gas monitor and connect gas measuring line (CO<sub>2</sub> line) to the AnaConDa and the gas monitor.
2. Connect the AnaConDa agent line with the red check valve to the AnaConDa syringe inside the syringe pump.
3. Remove the red protective cap from the AnaConDa and connect AnaConDa between patient and respirator Y-piece.
4. Program and apply a bolus of 1.5 ml for filling the AnaConDa agent supply line and the evaporator. Caution when administering bolus (do not overdose by using oversized bolus).

*Syringe pumps with programmable bolus rates are preferred so that venting the AnaConDa agent supply line won't result in overdosing. When using non-programmable syringe pumps be sure to administer an accurate bolus dose.*

### TIP

- Always use room temperature volatile anaesthetics.
- Avoid unnecessary manipulation (e.g. shaking) of bottle with volatile anaesthetic (pressurized bottle).
- Avoid heat sources which impact the stability of volatile anaesthetics (e.g. evaporation inside syringe).

## 4. Instructions for Use

### Filling and placing the syringe

1. Only use fill adapters by Sedana Medical with Isoflurane- or Sevoflurane bottles (standard thread or Sevoflurane Quickfill closure by Abbott).

The respective fill adapter is screwed to the bottle and may remain there until the contents have been used up (*the fill adapter for Quickfill bottles by Abbott must remain on the bottle until the bottle is completely empty / this avoids leaks*).

A built-in check valve inside the fill adapter prevents the volatile anaesthetic from escaping the bottle.

2. Remove syringe from AnaConDa packaging. Fill syringe with some air (approx. 20 ml).

3. Connect syringe to the respective fill adapter by firmly pressing and turning, and slowly and carefully fill syringe with liquid volatile anaesthetic by moving the syringe plunger back and forth a few times (do not generate too much excess/negative pressure).

After filling the syringe disconnect the syringe from the bottle, remove any air bubbles from the syringe. Close the syringe with the attached red cap.

4. Place the syringe into the syringe pump.

1



2



3



4



## Initialisation / AnaConDa set-up

1. Switch on the prepared gas monitor and connect gas measuring line (CO<sub>2</sub> line) to the AnaConDa and the gas monitor.



2. Connect the AnaConDa agent line with the red check valve to the AnaConDa syringe inside the syringe pump.



3. Remove the red protective cap from the AnaConDa and connect AnaConDa between patient and respirator Y-piece.



4. Program and apply a bolus of 1.5 ml for filling the AnaConDa agent supply line and the evaporator. Caution when administering bolus (do not overdose by using oversized bolus).


Start syringe pump at a rate of:

- 3 ml/h Isofluran or
  - 5 ml/h Sevofluran
- and adjust dose as necessary.



## Dosage calculation for end-tidal concentration (Fet)

The following dosage nomogram provides a reference value applicable to an end-tidal concentration of (Fet) 0.5Vol% Isoflurane or Sevoflurane.



Curve 1 represents the values during the first hour

Curve 2 represents the values starting with the 2nd hour

Initial syringe pump rates of Isoflurane and Sevoflurane have been established as follows based on studies and expert opinions:

**Isofluran: 3 ml/h**

**Sevofluran: 5 ml/h**

Sedation must be individually evaluated with each patient. Once the intended sedation level (sedation score) has been reached another adjustment should be made within the first hour of use by administering a bolus (0.5 ml) or modifying the syringe pump rate in accordance with additional or reduced sedation.

Please note the different MAC values of the various anaesthetic gases (see chapter volatile anaesthetics).

## Practical application

### Syringe changes

AnaConDa is a 24 hour single-patient device. The syringe may possibly (syringe empty) need to be replaced multiple times in a 24 hour period. Only use AnaConDa syringes.

1. Fill the new syringe with the volatile anaesthetic as instructed
2. Stop syringe pump
3. Disconnect AnaConDa agent line from the syringe in the syringe pump
4. Remove the empty syringe from the syringe pump
5. Place the new syringe into the syringe pump
6. Connect the AnaConDa agent line to the syringe
7. Restart syringe pump

### AnaConDa System change

AnaConDa is a 24 hour single-patient medical device. The system is changed every 24 hours.

### Suctioning

- Use a closed suction system.
- Activate the extraction key (O<sub>2</sub> oxygenation) on ventilator; wait for ventilator to pause → disconnect and suction. Use a swivel adapter with bronchoscopy cap.

### Using medication nebulisers

- Medication nebulisers are always attached before the AnaConDa, close to the patient.
- Ultrasound nebuliser do not impact the end-tidal gas concentration (F<sub>et</sub>).
- Jet air nebulisers require a pressurized air flow of up to 7 litres, which could result in a decrease of the end-tidal concentration (F<sub>et</sub>). Adjust sedation for this period (by increasing the syringe pump rate).

## Disconnection / connection

### **Always stop the perfusor when disconnecting!**

AnaConDa only releases anaesthetic gas when the flow is directed towards the patient.

- **Always disconnect** the AnaConDa at the Y-piece first, then from the patient
- **Always first connect** the AnaConDa to the patient first, then to the Y-piece

## Patient transport with AnaConDa

The reflexion properties of the AnaConDa allow the patient-side gas concentration to be maintained for some time before a drop in the concentration is noticeable.

Patients can be transported with AnaConDa. Sedation will be maintained at about the same level for approx. 25-35 minutes without continuing the infusion with the syringe pump. Based on experience a small bolus is administered shortly before patient transport. (Caution – avoid overdosing from bolus administration)

The syringe pump may also be included during transport.

## Temporary concentration increase

In some cases it may be necessary to quickly increase the sedation depth on a patient. In this case the bolus function of the syringe pump may be used (Caution – avoid overdosing from bolus administration). Bolus with approx. 0.2-0.3 ml Isoflurane or 0.4-0.6 ml Sevoflurane is sufficient to raise the  $F_{et}$  value by 0.2-0.4Vol%.

## 5. Warnings, Notices, Tips And Tricks

### Warnings

- Only use room-temperature anaesthetic gas (Isoflurane or Sevoflurane)!
- Avoid overdosing by incorrect bolus administration. Programmable syringe pumps are preferred!
- Never manually fill the agent line. Always use a syringe pump!
- Never close the AnaConDa with caps or other items whilst not connected to the patient!
- Always stop the syringe pump when disconnecting the AnaConDa!
- AnaConDa is a 24 hour single-patient product. Single use only!

### Notices

- Set the syringe setting on Becton Dickinson Plastipak or Sherwood Mo-noject 50ml on the syringe pump
- Set the syringe pump to maximum alarm tolerance (switch-off pressure to maximum)
- If necessary, remove any air or gas bubbles from the syringe.
- Always position AnaConDa with the black side facing up.
- Always keep the fill adapter for Sevoflurane Quickfill bottles connected to the bottle until the bottle is completely emptied. Repeated screwing on and unscrewing will result in bottle leaks.
- Please note the AnaConDa creates an additional dead space of approx. 100ml.
- Adjust the ventilation parameters if the CO<sub>2</sub> value rises too high.
- Always connect the medication nebulizer on the patient's side, before the AnaConDa
- Modifying the flow (ventilation parameters, compressed air nebulizer) impacts the end-tidal gas concentration.
- Always first disconnect the AnaConDa from the respirator (Y-piece).
- There may be various causes for incorrect or fluctuating gas monitor values (self-respiration, moisture in gas measuring line, AnaConDa position changes e.g. patient in Rotorest bed)
- AnaConDa has a very efficient HME (heat and moisture exchange) function. If necessary, remove any condensation on the patient end.

## Tips and tricks

### **Can AnaConDa be used for neonatal or paediatric purposes?**

The dead space of 100 ml limits the use of AnaConDa, we therefore recommend only using it on patients with a tidal volume of >350 ml. Use below 350 ml has been investigated. In this case the AnaConDa is built into the inspiratory tube and used purely as a vaporizer. The storage/recycle medium function is not applied. In small tidal volumes the consumption of volatile anaesthetic is not very high.

### **What happens if the black AnaConDa side is positioned upside-down/downward?**

Condensation could be drawn into the gas measuring line which will falsify the gas concentration measurement. Condensation on the hydrophobic bacteria/virus filter membrane may prevent air from passing. In addition the gas measurement values are false, as the heavy gas "drips" into the line instead of being drawn in showing falsely elevated levels.

### **Which gas measurement line should be used?**

Always use the measurement line recommended for the gas monitor. Incorrect lengths and types may impact the measurement. A special Nafion line (semi-permeable layer) can reduce the accumulation of condensation inside the water trap.

### **Why are end-tidal values ( $F_{et}$ ) higher than inspiratory values ( $F_i$ ) when using the AnaConDa?**

The AnaConDa performance differs from constant gas flow in an anaesthesia machine. The highly efficient storage medium (anaesthesia reflector) releases a large amount of anaesthetic gas to the patient at the start of inspiration, which would be shown as a so-called peak in the measured curve. The gas monitor, however, measures at the end of inspiration, when the majority of anaesthetic gas has already been fed to the patient.

This results in a lower than actual  $F_i$  value. The  $F_{et}$  value is indicated correctly.

### **High fluctuation in $CO_2$ values?**

When using a bypass flow monitor, check for water in the line. Check if the water trap is full. Remove any condensation.

## 6. Literature

### AnaConDa in intensive care units

1. Meiser A, Bellgardt M, Vogelsang H, Sirtl C, Weber T. Funktionsweise „des Anaesthetic Conserving Device“. Der Anaesthesist DOI 10.1007/s00101-010-1779-6
2. Soukup J, Schörrf K, Kubosch, Pohl C, Bomplitz M, Kompardt J. Sedation concepts with volatile anesthetics in the critically ill patients. *Journal of Critical Care* (2009) 24, 535-544.
3. Veismann U, Rehberg S, Rehberg J, Westphal M. Inhalationsanästhesie mittels „Anaesthetic Conserving Device“ zur Langzeitsedierung eines schwer sedierbaren Patienten.
4. Röhm K, Piper S, Boldt J. Prolonged Inhalational Sedation Using Sevoflurane: Evaluation of Inorganic Fluoride levels and Kidney Functions. *Adv. anesth crit care* Volume 1, Issue 2. 2009
5. Rohm KD, Mengistu A, Boldt J, Mayer J, Beck G, Piper SN. Renal integrity in Sevoflurane with the Anesthetic Conserving Device in the Intensive Care Unit – A comparison to intravenous propofol sedation. *Anesth Analg* 2009;108:1848-54.
6. Sackey PV, Martling CR, Carlswald C, Sundin O, Radell PJ. Short- and long-term follow-up of intensive care unit patients after sedation with isoflurane and midazolam – a pilot study. *Crit Care Med* 2008;36:801-6.
7. Rohm KD, Wolf M, Schollhorn T, Schellhaass A, Boldt J, Piper SN. Short-term sevoflurane sedation using the Anaesthetic Conserving Device after cardiothoracic surgery. *Intensive Care Med* 2008;34:1683-9.
8. Kompardt J, Scharff K, Kubosch K, Pohl C, Bomplitz M, Soukup J. Sedierung mit volatile Anästhetika auf der Intensivstation : Praktische Anwendung und derzeitige Erfahrungen mit dem AnaConDa-System. *Anaesthesist* 2008;57:1201-9.
9. Jung C, Granados M, Marsol P, Murat I, Gall O. Use of sevoflurane sedation by the AnaConDa device as an adjunct to extubation in a pediatric burn patient. *Burns* 2008;34:136-8.
10. Belda JF, Soro M, Badenes R, Meiser A, Garcia ML, Aguilar G, Marti FJ. The predictive performance of a pharmacokinetic model for manually adjusted infusion of liquid sevoflurane for use with the Anaesthetic-Conserving Device (AnaConDa): a clinical study. *Anesth Analg* 2008;106:1207-14.
11. Soukup J, Scharff K, Kubosch K, Rudzki M, Zielonka M, Radke J. Sedierung mit volatile Anästhetika auf der Intensivstation: Technische Umsetzung und derzeitige Erfahrungen. *Intensiv- und Notfallbehandlung* 2007;32:29-36.
12. Nickell EA, Benken I, Bartels U, Voelckel WG, Quintel M. AnaConDa als Ultima-Ratio-Therapie. Fallbericht einer chronisch obstruktiven Lungenerkrankung. *Anaesthesist* 2007;56:587-91.
13. Elliot S, Berridge JC, Mallick A. Use of the AnaConDa anaesthetic delivery system in ICU. *Anaesthesia* 2007;62:752-3.
14. Berton J, Sargentini C, Nguyen JL, Belii A, Beydon L. AnaConDa reflection filter: bench and patient evaluation of safety and volatile anesthetic conservation. *Anesth Analg* 2007;104:130-4.
15. Sackey PV, Martling CR, Radell PJ. Three cases of PICU sedation with isoflurane delivered by the 'AnaConDa'. *Paediatr Anaesth* 2005;15:879-85.
16. Sackey PV, Martling CR, Nise G, Radell PJ. Ambient isoflurane pollution and isoflurane consumption during intensive care unit sedation with the Anesthetic Conserving Device. *Crit Care Med* 2005;33:585-90.
17. Hanafy M. Clinical evaluation of inhalation sedation following coronary artery bypass grafting. *Eg J Anaesth* 2005;21:237-42.

18. Sackey PV, Martling CR, Granath F, Radell PJ. Prolonged isoflurane sedation of intensive care unit patients with the Anesthetic Conserving Device. *Crit Care Med* 2004;32:2241-6.

## **Volatile anaesthetics in intensive care units**

19. Meiser A, Laubenthal H. Inhalational anaesthetics in the ICU: theory and practice of inhalational sedation in the ICU, economics, risk-benefit. *Best Pract Res Clin Anaesthesiol* 2005;19:523-38.
20. Meiser A, Sirtl C, Bellgardt M, Lohmann S, Garthoff A, Kaiser J, Hugler P, Laubenthal HJ. Desflurane compared with propofol for postoperative sedation in the intensive care unit. *Br J Anaesth* 2003;90:273-80.
21. Tanigami H, Yahagi N, Kumon K, Watanabe Y, Haruna M, Matsui J, Hayashi H. Long-term sedation with isoflurane in postoperative intensive care in cardiac surgery. *Artif Organs* 1997;21:21-3.
22. Osborne MA, Eddleston JM, McNicol W. Inorganic fluoride concentration after long-term sedation with isoflurane. *Intensive Care Med* 1996;22:677-82.
23. Kong KL. Inhalational anaesthetics in the intensive care unit. *Crit Care Clin* 1995;11:887-902.
24. Spencer EM, Willatts SM. Isoflurane for prolonged sedation in the intensive care unit; efficacy and safety. *Intensive Care Med* 1992;18:415-21.
25. Millane TA, Bennett ED, Grounds RM. Isoflurane and propofol for long-term sedation in the intensive care unit. A crossover study. *Anaesthesia* 1992;47:768-74.
26. Johnston RG, Noseworthy TW, Friesen EG, Yule HA, Shustack A. Isoflurane therapy for status asthmaticus in children and adults. *Chest* 1990;97:698-701.
27. Kong KL, Willatts SM, Prys-Roberts C. Isoflurane compared with midazolam for sedation in the intensive care unit. *Brmj* 1989;298:1277-80.
28. Kofke WA, Young RS, Davis P, Woelfel SK, Gray L, Johnson D, Gelb A, Meeke R, Warner DS, Pearson KS, et al. Isoflurane for refractory status epilepticus: a clinical series. *Anesthesiology* 1989;71:653-9.

## **Sedation in general**

29. Kessler P, Martin J. Analgosedierung. Umsetzung der S2e-Leitlinien dient der Prozessoptimierung. *Anesthesiol Intensivmed Notfallmed Schmerzther* 2008;43:38-43.
30. Kress JP, Vinayak AG, Levitt J, Schweickert WD, Gehlbach BK, Zimmerman F, Pohlman AS, Hall JB. Daily sedative interruption in mechanically ventilated patients at risk for coronary artery disease. *Crit Care Med* 2007;35:365-71.
31. Martin J, Franck M, Fischer M, Spies C. Sedation and analgesia in German intensive care units: how is it done in reality? Results of a patient-based survey of analgesia and sedation. *Intensive Care Med* 2006;32:1137-42.
32. Kress JP, Hall JB. Sedation in the mechanically ventilated patient. *Crit Care Med* 2006;34:2541-6.
33. Martin J, Parsch A, Franck M, Wernecke KD, Fischer M, Spies C. Practice of sedation and analgesia in German intensive care units: results of a national survey. *Crit Care* 2005;9:R117-23.
34. Schaffrath E, Kühlen R, Tonner PH. Analgesie und Sedierung in der Intensivmedizin. *Anaesthesist* 2004;53:1111-30.

## 7. Sedation Scoring Systems

Scoring systems helps the intensive care personnel to quite accurately assess the patient's sedation level through the use of clinical criteria. The most frequently used sedation scores are:

### Ramsey Sedation Scale

Score	Description	Assessment
0	Awake, oriented	Awake
1	Agitated, restless, anxious	Too flat
2	Awake, cooperative, respiratory tolerance	Adequate
3	Asleep but cooperative (opens eyes in response to loud verbal commands or touch)	Adequate
4	Deep sedation (doesn't open eyes in response to loud verbal commands or touch, but promptly reacts to pain stimulation)	Adequate
5	Anaesthesia (lethargic reaction to pain stimulation)	Deep
6	Deep coma (no reaction to pain stimulation)	Too deep

## The Richmond Agitation Sedation Scale (RASS)

Score	Term	Description
+4	Combative	Combative or aggressive, immediate danger to staff
+3	Very agitated	Pulls on or removes tubes, catheter, etc., or aggressive toward staff
+2	Agitated	Frequent non-purposeful movement or patient-ventilator dyssynchrony
+1	Restless	Anxious but movements are not aggressive or vigorous
0	Alert and calm	
-1	Drowsy	Not fully alert but has sustained awakening (more than 10 seconds), eye contact to voice
-2	Light sedation	Briefly awakens (less than 10 seconds) with eye contact to voice
-3	Moderate sedation	Movement in response to voice, but no eye contact
-4	Deep sedation	No response to voice, but movement in response to physical stimulation
-5	Unarousable	No response to voice or physical stimulation

## Material List

Product/component	Material	Residue after burning
-------------------	----------	-----------------------

### AnaConDa System

Housing	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O
Evaporator	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O
Carbon filter	Carbon	CO <sub>2</sub>
Bacteria and virus filter	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O
Agent line	Polyethylene	CO <sub>2</sub> , H <sub>2</sub> O
Line connector	Polyethylene	CO <sub>2</sub> , H <sub>2</sub> O
Hot glue	Olefin co-polymer	CO <sub>2</sub> , H <sub>2</sub> O
Luer cap	Polyethylene	CO <sub>2</sub> , H <sub>2</sub> O
Locking plug	Polycarbonate	CO <sub>2</sub> , H <sub>2</sub> O
Label	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O

### Syringe

Barrel	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O
Plunger	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O
Piston	Rubber	CO <sub>2</sub> , H <sub>2</sub> O
Adhesive	Acrylic U.V.	CO <sub>2</sub> , H <sub>2</sub> O
Lubricant	Silicone	SiO <sub>2</sub>
Label	Polypropylene	CO <sub>2</sub> , H <sub>2</sub> O

### Packaging

Over layer	Paper	CO <sub>2</sub> , H <sub>2</sub> O
Tray	Polyethylene terephthalate	CO <sub>2</sub> , H <sub>2</sub> O

## Operating Conditions

Anaesthetic agents:	Isoflurane or Sevoflurane
Tidal volume operating range:	minimum 350 ml
Resistance @ 60 l/min:	2.5 cm H <sub>2</sub> O (250 Pa)
Moisture loss at 0.75 l x 12 breaths/min:	5 mg/l
Humidification output (calc.):	30 mg/l
Moisture loss at 1.0 l x 10 breaths/min:	7 mg/l
Humidification output (calc.):	29 mg/l
Bacterial filtration efficiency:	99.999%
Viral filtration efficiency:	99.98%
Dead space:	approx. 100 ml
Weight:	50 g
Connectors per ISO 5356:	15F/22M - 15 M
Gas monitoring port:	Female Luer lock
Agent line length:	2200 mm

## Ordering Information

Item no.	Description
26000	AnaConDa System with syringe
26022	AnaConDa syringe, single
26064	Filling adapter 26064 Isoflurane and Sevoflurane with standard threads
26042	Filling adapter 26042 for Sevoflurane Quickfill by Abbott
26072	Accessory set for rest gas filter
ZE0000050	Contrafluran anaesthetic gas filter
ZE0000051	Contrafluran standard mount
ZE0000052	Contrafluran Sensoflurane mount

## Contributors

### **Mr. OA Dr. med. Andreas Meiser**

Saarland University Hospital  
Clinic for Anaesthesiology, Intensive Care and Pain Management  
Kirrberger Straße 100  
D-66421 Homburg/Saar

### **Mrs. OÄ Priv.-Doz. Dr. med. Kerstin D. Röhm**

Clinical Centre of the City of Ludwigshafen gGmbH  
Clinic for Anaesthesiology and Surgical Intensive Care  
Bremsersstraße 79  
67063 Ludwigshafen

### **Mr. OA Dr. med. Jens Soukup**

University Hospital Halle (Saale)  
University Clinic for Anaesthesiology and Surgical Intensive Care  
Ernst-Grube-Str. 40  
06120 Halle (Saale)

## Legal Basis

The AnaConDa<sup>®</sup> system is a CE-certified medical device. Use of the system and its components during “inhalation sedation” in intensive care patients is based on the Medical Devices Act (MPG, Medizinproduktegesetz in its version published August 7, 2002 (BGBl. I pg. 3146), last modified with Article 1 of the Act dated June 14, 2007 (BGBl. I pg. 1066).

When using Sevoflurane or Isoflurane please follow the respective manufacturer expert information.



# SEDANAMEDICAL

the AnaConDa® technology people



3 000 025-1009

**Address** Kungsgatan 62, SE-753 18 Uppsala, Sweden  
**Phone** +46-(0)18-18 66 27, Fax +46-(0)18-10 15 80  
**Email** [info@sedanamedical.com](mailto:info@sedanamedical.com)  
**Internet** [www.sedanamedical.com](http://www.sedanamedical.com)