

Say hello to the all new AnaConDa®. Sedation has never been this easy.

A new valve on the AnaConDa side prevents leakage when not connected to the syringe

Now with unique connectors between the AnaConDa and the syringe to eliminate the risk for connecting anything with a luer connector

The new filling adaptor comes with a new dimension to fit the new non-luer connector



.....► You can download the latest Instruction for Use from www.sedanamedical.com

SEDANAMEDICAL
the AnaConDa® technology people

The all new AnaConDa[®]

Anaesthetic Conserving Device



AnaConDa is a disposable administration system for isoflurane and sevoflurane. It opens up the possibility to use volatile agents in an easy and safe way without the need for an anaesthesia machine. It can be used for sedation in the ICU and for anaesthesia.

Main advantages:

Short and predictable wake up time

- extubation can easily be planned

No development of tolerance

- the same dose can normally be used throughout the treatment
- low concentrations usually needed ~ 0,3% (isoflurane)

More reliable route of administration and elimination

- Not depending on hepatic metabolism and/or renal excretion for elimination

No withdrawal symptoms

Sedation levels are easy to control

- Alveolar concentration can be continuously monitored with a gas monitor.

“Prolonged isoflurane sedation of intensive care unit Patients with the Anaesthetic Conserving Device”

Sackey P, Martling C-R, et al. Crit Care Med 2004 Vol. 32, No.11

“The Anaesthetic Conserving Device allows titratable administration of isoflurane without costly equipment and can safely be managed by the nursing staff”

Sackey et al: Crit.Care.Med 2004; 32: 2241-2246

“It has been shown that using isoflurane is an environmentally safe method of sedation when used together with the AnaConDa”

Sackey et al: Crit Care Med 2005; Vol 33 585-590

	Isoflurane	Midazolam
Time to Extubation (mins)	10±5	250±270
Time to follow verbal command (mins)	10±8	110±130

SEDANA MEDICAL

the AnaConDa[®] technology people

Go to our web site www.sedanamedical.com to find your local dealer and to download all material you need to be able to use AnaConDa[®]